

## T-Link Cell Phone Coupling Device



This device works with your hearing aid telecoil to help with hearing cell phones.

### Features:

- Gives you hands-free operation
- Built in microphone
- Helps avoid interference from your cell phone and can improve clarity of sound
- Lets you hear in both ears so you can better understand a conversation
- Its standard 3/32 inch (2.5 mm) plug can be used with most cell phones.
- Does not require any batteries
- Size is 50" – 34" from the plug to the "Y" section of the cord
- For a longer cord, Radio Shack has a 6-foot extension